

CDI020190786

學與教策略系列：運用資訊科技工具以進行科學探究學習活動

**Simple Colorimeter constructed by Micro:bit and Smartphone
(ChemEye App)**


I. Components of Device



Micro:bit



Solderless Breadboard



RGB-LED light




Electric Wire



Extension Board




3D-Printing Box



Resistors (220-Ohm)

Set-up:




| | | | |
|---------------|-------------------------------|----------------------------------|--|
| RGB-LED Light | Color of Wire used to connect | Circuit of Solderless Breadboard | Pin of Extension board Chosen to connect |
| R | Red | | Pin0 |
| GND | Brown | | / |
| G | Green | | Pin1 |
| B | Blue | | Pin2 |

- II. Prepare the following concentrations of Standard Solution (Tartrazine 檸檬黃).
0ppm(Blank), 20ppm, 40ppm, 60ppm, 80ppm, 100ppm
- III. Download ChemEye App from Apple Store(IOS) or Google Play(Android)



- IV. Select the most suitable RGB value for Tartrazine.



- V. Coding (<https://www.microbit.org/code/>)

| Color | Value | Emitted Color | Absorbed Color |
|-------|-------|---------------|----------------|
| R | 255 | — | — |
| G | 255 | — | — |
| B | 255 | — | — |

forever

analog write pin P0 to map 59 from low 0 high 255 to low 0 high 1023

analog write pin P1 to map 86 from low 0 high 255 to low 0 high 1023

analog write pin P2 to map 239 from low 0 high 255 to low 0 high 1023

VI. Different concentrations of standard solution are transferred into the cuvette and are measured by ChemEye App. The readings can be collected.

| | Blank | 20ppm | 40ppm | 60ppm | 80ppm | 100ppm | Unknown |
|---|-------|-------|-------|-------|-------|--------|---------|
| R | | | | | | | |
| G | | | | | | | |
| B | | | | | | | |

Choose the data of _____ (R / G / B) value for plotting calibration curve.

| | |
|--------------------------------------|---------------------------|
| Transmittance: $T = \frac{P_n}{P_0}$ | Absorbance: $A = -\log T$ |
|--------------------------------------|---------------------------|

| | Blank | 20ppm | 40ppm | 60ppm | 80ppm | 100ppm | Unknown |
|---------------|-------|-------|-------|-------|-------|--------|---------|
| Transmittance | | | | | | | |
| Absorbance | | | | | | | |

Using the given graph paper to plot the calibration curve and to **Find out the concentration of Unknown solution.**

VII. The concentration of unknow solution = _____.